

# Domilia Sit-to-Stand

## LIA3000



### Technical sheet

Product code	LIA3000
Dimensions	- Base : Length: 24" ½ x Width: 23 ½" - Upper support bar : Height: 48" Width: 19" - Profondeur : 15"
Materials	- Base and upper support bar : Steel - Leg support : Foam
Other characteristics	- 5" Casters with integrated wheels - Adjustable height and width for leg support - Non-slip base
Product weight	20.4 kg (45 lbs), Base 11 kg (25 lbs), Lateral bar: 9 kg (20 lbs)
Assembly required	Yes – 1 simple step
Warranty	Manufacturer warranty – 1 year
Maximum weight capacity	400 lbs



#### Recommended use :

##### Sit-to-Stand aid only :

- Someone with upper body muscle tone
- Ability to stand up, but difficulty taking a few steps
- Good standing balance

##### Sit-to-Stand aid with belt

- Someone with upper body muscle tone
- Ability to stand up, but difficulty taking a few steps
- Weight bearing on legs, lack of stamina when standing
- Fall risks

#### Environments for product use :

- Home care
- Long-term care home and facility
- Hospitals
- Rehabilitation center



Innovate - Prevent - Communicate